

Title:	Replogle Tube Insertion and Care		
Authored by:	London Neonatal Transfer Service	Reviewed date:	December 2014
Reviewed by:	London Neonatal Transfer Service	Next review date:	December 2016

Aim is to maintain upper airway and keep the upper pouch free of secretions.

Equipment

- Replogle tube size 8 or 10 French
- Size 10 Fr nasogastric / orogastric tube
- Low pressure wall suction unit
- Suction tubing
- 5mls syringe with occlusion cap (red cap)
- 2 ampoules of 0.9% saline
- 5 and 10ml syringes
- Scissors
- Duoderm, Tegaderm or any clear transparent dressing
- Monitoring equipment: ECG, respirations and saturations
- Spare suction catheters and Replogle tubes

The baby should be fully assessed and monitoring of vital signs implemented.

Procedure

- Explain procedure to parents to alleviate anxiety
- Wash hands and put on non-sterile gloves to prevent cross infection
- Suction nasal passages and oropharynx if required to clear airway.
- Pass Replogle tube via nostril until resistance is felt, then withdraw slightly. Fix tube to face using Duoderm and Tegaderm or any clear transparent dressing, ensuring as little distortion of nostril as possible. Ensure tube is well fixed and minimise trauma to nostril and mucosa of pouch. Record the tube size and length at nose in the notes.
- Tube should then be attached to low level wall suction at 5-10 kPa. High level suction can traumatise the oesophageal pouch and impede surgical repair.
- If secretions are thick, it may be necessary to increase frequency of flushes. Do not increase the suction pressure.
- Instil 0.5–1.0ml of 0.9% saline every 15 min to maintain patency. Saline should flow back quickly, slow flowback indicates the tube may be blocking. Increase the frequency of flushing and / or reposition the tube.
- Observe baby for signs of blockage or displacement of tube i.e.: flushes not being returned, coughing, cyanosis and desaturation. Baby should be observed for signs of distress i.e.: tachypnoea, recessions and nasal flaring.
- Spare Replogle tube of correct size, tape measure and scissors to be kept close to the baby.
- Suction bag and tubing needs to be changed following the transfer.

Note: If continuous low pressure suction is unavailable, then the other option is to pass wide bore naso/oro gastric tube and then suction regularly with a syringe.